




FILLMED
LABORATOIRES

ART FILLER®
VOLUME

USE

ART FILLER®
VOLUME

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DESCRIPTION

VOLUME is a viscoelastic cross-linked hyaluronic acid gel of non-animal origin gradually absorbed over time. It is colourless, transparent, sterile, non-pyrogenic and saline and contains a 0.3 % volume of lidocaine hydrochloride for its anaesthetic properties. It is available in a pre-filled and single-use 1.2 mL syringe. Each pack contains 2 VOLUME syringes, 4 sterile and single-use 27G½" needles intended for injection of VOLUME, a product leaflet and 4 traceability labels. In order to guarantee traceability of the product, one of the two labels must be affixed to the patient file and the other must be transmitted to the patient.

INGREDIENTS

Cross-linked hyaluronic acid	25 mg
Lidocaine hydrochloride	3 mg
Phosphate buffer solution, pH 7.2	QS 1 g
A syringe contains 1.2 mL of VOLUME gel.	

STERILISATION

The content of the VOLUME syringes is sterilised using moist heat. The 27G½" needles are sterilised using radiation.

INTENDED USE / INDICATIONS

ART FILLER® products are intended to restore changes to skin structure caused by ageing: fills in lines and creases in the facial skin, restores volume. VOLUME filler is an injectable implant indicated to restore volume to the mid-face, jawline and temples, via subcutaneous, supraperiosteal or deep dermis injections. The inclusion of lidocaine aims to reduce painful sensations for the patient during treatment.

PERFORMANCE

In a multicentre prospective study with ART FILLER® VOLUME for medio-facial restoration, 66 % of the patients maintained a clinically significant improvement 12 months after treatment, and 37 %, 18 months after treatment.

CONTRAINDICATIONS

VOLUME must not be injected:

- For correction of superficial and fine lines.
- Into the eye area (eyelids, crows' feet, under eye area), the glabellar area or into the lips.
- Into blood vessels.
- Into muscles.
- Into an area where a non-absorbable filler implant has already been injected.
- Do not over-correct.

VOLUME must not be injected:

- In patients with known hypersensitivity to hyaluronic acid, lidocaine and amide type local anesthetics.
- In patients with a medical history of auto-immune illness or receiving immunotherapy.
- In patients suffering from epilepsy that is not controlled with treatment.
- In patients affected by porphyria.
- In patients with a tendency to develop hypertrophic scars.
- In patients with medical history of recurring sore throats related to rheumatic fever localised in the heart.
- In pregnant or breastfeeding women.
- In children.
- In areas that have inflamed and/or contagious skin lesions (acne, herpes, etc.).
- Immediately after or before laser treatment, deep chemical peeling or dermabrasion.

**WARNINGS**

- VOLUME is not indicated for injections other than subcutaneous, supraperiosteal or deep dermis injections. The technique and the depth of injection vary according to the treatment area.
- Do not inject into blood vessels, bones, tendons, ligaments, muscles or beauty spots.
- Unintentional injection of soft tissue fillers into blood vessels in the face can result in rare, but serious side effects such as embolization which can cause vision impairment, blindness, stroke and damage and/or skin and underlying facial structures necrosis. These rare cases of blood vessels embolization are mostly reported in glabella, in and around the nose, forehead, and periorbital region.
- Check the expiry date and check that the packaging is intact before use. Do not use a syringe with opened or shifted tip cap within the blister.
- Do not reuse. The re-use of a product bears a risk (e.g. cross-contamination) for the patient.
- Do not re-sterilise.
- Dispose of the syringe and the remaining product after use. The needle must be disposed of in a container provided for this purpose. Refer to the guidelines in force for their disposal.
- Never try to straighten a bent needle. Dispose of it and replace it.

PRECAUTIONS PRIOR TO USE

- The use of VOLUME is reserved for doctors trained in injection techniques.
- There are no clinical data available in terms of efficacy and tolerance of the VOLUME injections for an area that has already been treated with another filler product.
- There are no clinical data available in terms of efficacy and tolerance of the VOLUME injections for Fitzpatrick phototypes V and VI.
- There are no clinical data available in terms of efficacy and tolerance of the VOLUME injections in patients who have a medical history of severe multiple allergies or anaphylactic shock. Therefore, the doctor must decide about the indication on a case by case basis according to the nature of the allergy

and (s)he must specifically monitor patients who present a risk. In particular, (s)he may decide to offer a double test or adapted preventative treatment prior to any injection.

- VOLUME must be used with precaution in patients with cardiac conduction conditions.
- VOLUME must be used with a great deal of precaution in patients suffering from hepatocellular insufficiency with coagulation disorders as well as patients receiving treatment with medicinal products that reduce or inhibit hepatic metabolism which are prone to result in coagulation disorders.
- The injection volume depends on the correction required and is under the discretion of the practitioner. The recommended dose per treatment site is 1.2 mL. A maximum of two treatment sessions per year is recommended. Possible touch-up sessions enable the sought-after degree of correction to be maintained. The recommended maximal dose per year is 15 mL. There is no data available regarding the safety of injecting greater amount than 15 mL of ART FILLER® dermal fillers per 70 kg body mass per year.
- Patients must receive the following advices:
 - Avoid taking aspirin and vitamin C and/or vitamin E in high doses the week prior to the injection.
 - Patient receiving anticoagulant treatment must be warned of the increased risk of bruising and bleeding during injection.
 - Do not apply any make-up at all for 12 hours after the injection.
 - Avoid exposure to extreme temperatures (extreme cold, sauna, steam room) as well as prolonged sun exposure or ultraviolet light for 2 weeks following the injection.
- If the needle is obstructed, do not increase the pressure on the tip of the plunger. Stop the injection and use a new needle.
- Doctors are reminded to bear in mind that this product contains lidocaine and they must take this into account.
- Sports participants are reminded that this product contains an active ingredient that may lead to a positive result in any drug tests carried out.

INCOMPATIBILITES

Hyaluronic acid is incompatible with quaternary ammonium compounds, such as benzalkonium chloride solutions. This is why VOLUME must never come into contact with medical and surgical instruments treated with this type of product.

ADVERSE EFFECTS

The doctor must inform the patient that potential adverse effects exist that are related to implantation of this device occurring immediately or with a delay. Among these (non-exhaustive list) :

- Inflammatory reactions (redness, oedema, erythema) that may be combined with itching, and soreness when pressed may occur after the injection. These reactions may last for up to a week. In a limited number of cases, 1,6 % (1/61) of the subjects in the clinical study, slight pain, spontaneous or at palpation can last more than a week.
- Bruising.
- Change of colour or discolouring of skin in the injection area.
- Hardening or nodules in the area of the injection.
- Poor filling effectiveness or poor filling effect.
- Cases of necrosis of the glabella, abscesses, granuloma, peripheral facial paralysis and immediate or delayed hypersensitivity have been reported after injections of hyaluronic acid and/or lidocaine. These potential risks should be taken into account.
- The patient must inform the doctor as soon as possible if an inflammatory reaction persists for longer than one week or if any other side effects occur. The doctor will treat these appropriately.
- The distributor and/or manufacturer must be alerted to any other adverse effect related to the injection of VOLUME.

INSTRUCTIONS FOR USE

VOLUME is intended to be injected slowly into the deep dermis or supraperiosteal or subcutaneously a doctor trained in injection techniques. The technical skill of the treatment is essential for its success. This device must be used by practitioners who have carried out specific training for restoring volume. A good knowledge of the anatomy and physiology of the area to be treated is necessary. Before beginning the treatment, the patient must be informed of the indications of the device, its contraindications, its incompatibilities and its potential adverse effects. The area to be treated must be thoroughly disinfected before the injection. Use the 27G½" needle provided with the syringe. Remove the stopper of the syringe by pulling it out as indicated in diagram 1. Insert the needle correctly and firmly into the tip of the pre-filled syringe, twist it on carefully and ensure that the needle is in the correct position according to diagrams 2 and 3. Remove the cap of the needle (diagram 4) and slowly inject it into the dermis using the appropriate injection technique. The quantity injected will depend on the area to be treated. After injection, the doctor may massage the area lightly in order to distribute the product evenly. Graduations printed on the syringe label are intended to help the practitioner and cannot be considered as a measurement device.

STORAGE CONDITIONS

Store between 2°C and 25°C, away from direct sunlight.

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LABORATOIRES FILL-MED MANUFACTURING S.A.